

# Analytical Standards and CRMs for Food and Beverage Analysis

**Fluka**  
Analytical



**SUPELCO**  
Solutions within.™

 **Cerilliant**  
Analytical Reference Standards

 **RTC**

**LABCO L.L.C.**

**SIGMA-ALDRICH**

## Serving all your needs for Food and Beverage Analysis

Sigma-Aldrich brings decades of experience in the production and supply of analytical standards and reference materials for food and beverage analysis.

By staying in close contact with regulatory bodies and our customers, we maintain a range of standards for new and updated analytical methodologies. This flyer gives an overview of our current product offering.

We add continuously new products to the portfolio, in order to keep up-to-date with new developments in food and beverage analysis. Please visit our website at [sigma-aldrich.com/standardsfood](https://sigma-aldrich.com/standardsfood) for more information.



### Standards and Reference Material Categories

**Analytical Standards** are sufficiently pure, stable and homogeneous to be used for the calibration of an apparatus or assessment of analytical method performance. We provide a certificate of analysis for all analytical standards. Fluka Analytical Standards generally conform to or exceed the ISO Definition of a Reference Material.

**Certified Reference Materials (CRMs)** meet the requirements of ISO Guides 30 – 35 and are accompanied by a certificate stating the traceable value and its uncertainty. As part of the ISO Guide 34 requirements, homogeneity and stability of the materials are extensively tested and contribute to the uncertainties.

**Our in-house produced TraceCERT® CRMs** use high performance quantitative NMR (HP-qNMR®) under ISO/IEC 17025 to establish direct traceability to NIST SRM. On our webpage at [sigma-aldrich.com/organiccrm](https://sigma-aldrich.com/organiccrm) you can find more information and the complete product portfolio.

#### Matrix Materials from IRMM

The Institute for Reference Materials and Measurements (IRMM) provides a large number of CRMs for food and agriculture. These products are either pure materials (or synthetic mixtures) or matrix materials. All products are delivered with a certificate, which carries a certified value with its uncertainty traceable either to a SI unit or an internationally accepted reference.

### Custom Standards

If the standard you need is not available from our stock products, or if you need special packaging or testing, we can prepare custom standards according to your specific requirements.

To request a quote for an **organic custom standard**, please contact us at [customstandards@sial.com](mailto:customstandards@sial.com) or complete an online custom standards quote request at [sigma-aldrich.com/standards](https://sigma-aldrich.com/standards) for an **inorganic custom standard** for ICP or IC, please use our convenient online custom standards platform at [sigma-aldrich.com/csp](https://sigma-aldrich.com/csp)

### How to Find the Right Standard?

#### Online Standards Catalogue

Browse our entire line starting from our standards starting page [sigma-aldrich.com/standards](https://sigma-aldrich.com/standards) (click: Analytical Standards Product Catalogue).

#### Standards Explorer.

This search engine [sigma-aldrich.com/standards\\_explorer](https://sigma-aldrich.com/standards_explorer) provides a fast and efficient access to our extensive product line using keyword searching with name, CAS number, molecular formula or official methods.

## Food and Beverage Standards

Our portfolio of analytical standards and certified reference materials (CRMs) for food and beverage analysis comprise the following product groups:

### Neats and Solutions

Amino Acids, Peptides and Proteins

Antibiotics

Antioxidants

Aroclors / PCBs / Dioxins

Biotoxins / Mycotoxins

Carbohydrates

Carotenoids

DNA Standards  
(Food Authenticity Testing)

FA/FAMEs

Flavors and Fragrances

Food Color Additives

Food Residuals

Hormones

Inorganic Standards for AAS and ICP

Ion Standard Solutions

Lipids, Oils, Sterols

Natural compounds

NODGE/BADGE

Nutraceuticals

Organic Acids

Organic Neat CRMs

Pesticides and Metabolites

Phytopharma

Preservatives

Sweeteners

Veterinary Drugs and Metabolites

Vitamins

Water Standards for Karl Fischer Titration

### TraceCERT CRMs

Amino Acids

Antibiotics

FA/FAME

Mycotoxins

Organic Acids

Pesticides

Polyaromatic Hydrocarbons (PAHs)

### Matrix Materials

GMO Reference Materials  
(Genetically modified Organisms)

Food and Agriculture CRMs

### Biological Reference Materials

Microorganisms (CFU)

### Physical Properties

Color Reference Solutions

Density

Melting Point

Particle Size

pH Calibration / Buffers Redox

Viscosity

Find our products at [sigma-aldrich.com/standardsfood](https://sigma-aldrich.com/standardsfood)

**Expert Technical Service** is only a phone, FAX or email away. If you need help choosing the right standard or reference material, our technical service representatives will be happy to help you with:

- Product selection
- Obtaining certificates of analysis online
- Using a standard or reference material
- Troubleshooting product performance

**EurTechServ@sial.com**  
techservice@sial.com (US and Canada)

## Sigma-Aldrich Worldwide Offices

### Argentina

Free Tel: 0810 888 7446  
Tel: (+54) 11 4556 1472  
Fax: (+54) 11 4552 1698

### Australia

Free Tel: 1800 800 097  
Free Fax: 1800 800 096  
Tel: (+61) 2 9841 0555  
Fax: (+61) 2 9841 0500

### Austria

Tel: (+43) 1 605 81 10  
Fax: (+43) 1 605 81 20

### Belgium

Free Tel: 0800 14747  
Free Fax: 0800 14745  
Tel: (+32) 3 899 13 01  
Fax: (+32) 3 899 13 11

### Brazil

Free Tel: 0800 701 7425  
Tel: (+55) 11 3732 3100  
Fax: (+55) 11 5522 9895

### Canada

Free Tel: 1800 565 1400  
Free Fax: 1800 265 3858  
Tel: (+1) 905 829 9500  
Fax: (+1) 905 829 9292

### Chile

Tel: (+56) 2 495 7395  
Fax: (+56) 2 495 7396

### China

Free Tel: 800 819 3336  
Tel: (+86) 21 6141 5566  
Fax: (+86) 21 6141 5567

### Czech Republic

Tel: (+420) 246 003 200  
Fax: (+420) 246 003 291

### Denmark

Tel: (+45) 43 56 59 00  
Fax: (+45) 43 56 59 05

### Finland

Tel: (+358) 9 350 9250  
Fax: (+358) 9 350 9255

### France

Free Tel: 0800 211 408  
Free Fax: 0800 031 052  
Tel: (+33) 474 82 28 88  
Fax: (+33) 474 95 68 08

### Germany

Free Tel: 0800 51 55 000  
Free Fax: 0800 64 90 000  
Tel: (+49) 89 6513 0  
Fax: (+49) 89 6513 1160

### Hungary

Ingyenes telefonszám: 06 80 355 355  
Ingyenes fax szám: 06 80 344 344  
Tel: (+36) 1 235 9063  
Fax: (+36) 1 269 6470

### India

#### Telephone

Bangalore: (+91) 80 6621 9400  
New Delhi: (+91) 11 4358 8000  
Mumbai: (+91) 22 2570 2364  
Hyderabad: (+91) 40 4015 5488  
Kolkata: (+91) 33 4013 8003

#### Fax

Bangalore: (+91) 80 6621 9650  
New Delhi: (+91) 11 4358 8001  
Mumbai: (+91) 22 4087 2364  
Hyderabad: (+91) 40 4015 5488  
Kolkata: (+91) 33 4013 8000

### Ireland

Free Tel: 1800 200 888  
Free Fax: 1800 600 222  
Tel: (+353) 402 20370  
Fax: (+353) 402 20375

### Israel

Free Tel: 1 800 70 2222  
Tel: (+972) 8 948 4100  
Fax: (+972) 8 948 4200

### Italy

Tel: (+39) 02 3341 7310  
Fax: (+39) 02 3801 0737

### Japan

Tel: (+81) 3 5796 7300  
Fax: (+81) 3 5796 7315

### Korea

Free Tel: (+82) 80 023 7111  
Free Fax: (+82) 80 023 8111  
Tel: (+82) 31 329 9000  
Fax: (+82) 31 329 9090

### Malaysia

Tel: (+60) 3 5635 3321  
Fax: (+60) 3 5635 4116

### Mexico

Free Tel: 01 800 007 5300  
Free Fax: 01 800 712 9920  
Tel: (+52) 722 276 1600  
Fax: (+52) 722 276 1601

### The Netherlands

Free Tel: 0800 022 9088  
Free Fax: 0800 022 9089  
Tel: (+31) 78 620 5411  
Fax: (+31) 78 620 5421

### New Zealand

Free Tel: 0800 936 666  
Free Fax: 0800 937 777  
Tel: (+61) 2 9841 0555  
Fax: (+61) 2 9841 0500

### Norway

Tel: (+47) 23 17 60 00  
Fax: (+47) 23 17 60 10

### Poland

Tel: (+48) 61 829 01 00  
Fax: (+48) 61 829 01 20

### Portugal

Free Tel: 800 202 180  
Free Fax: 800 202 178  
Tel: (+351) 21 924 2555  
Fax: (+351) 21 924 2610

### Russia

Tel: (+7) 495 621 5828  
Fax: (+7) 495 621 5923

### Singapore

Tel: (+65) 6779 1200  
Fax: (+65) 6779 1822

### Slovakia

Tel: (+421) 255 571 562  
Fax: (+421) 255 571 564

### South Africa

Free Tel: 0800 1100 75  
Free Fax: 0800 1100 79  
Tel: (+27) 11 979 1188  
Fax: (+27) 11 979 1119

### Spain

Free Tel: 900 101 376  
Free Fax: 900 102 028  
Tel: (+34) 91 661 99 77  
Fax: (+34) 91 661 96 42

### Sweden

Tel: (+46) 8 742 4200  
Fax: (+46) 8 742 4243

### Switzerland

Free Tel: 0800 80 00 80  
Free Fax: 0800 80 00 81  
Tel: (+41) 81 755 2828  
Fax: (+41) 81 755 2815

### United Kingdom

Free Tel: 0800 717 181  
Free Fax: 0800 378 785  
Tel: (+44) 1747 833 000  
Fax: (+44) 1747 833 313

### United States

Toll-Free: 800 325 3010  
Toll-Free Fax: 800 325 5052  
Tel: (+1) 314 771 5765  
Fax: (+1) 314 771 5757

### Vietnam

Tel: (+84) 3516 2810  
Fax: (+84) 6258 4238

### Internet

[sigma-aldrich.com](http://sigma-aldrich.com)

*Enabling Science to  
Improve the Quality of Life*

Order/Customer Service (800) 325-3010 • Fax (800) 325-5052  
Technical Service (800) 325-5832 • [sigma-aldrich.com/techservice](http://sigma-aldrich.com/techservice)  
Development/Custom Manufacturing Inquiries **SAFC**® (800) 244-1173  
Safety-related Information [sigma-aldrich.com/safetycenter](http://sigma-aldrich.com/safetycenter)

World Headquarters  
3050 Spruce St.  
St. Louis, MO 63103  
(314) 771-5765  
[sigma-aldrich.com](http://sigma-aldrich.com)

©2012 Sigma-Aldrich Co. All rights reserved. SIGMA, SAFC, SIGMA-ALDRICH, ALDRICH, FLUKA, and SUPELCO are trademarks belonging to Sigma-Aldrich Co. and its affiliate Sigma-Aldrich Biotechnology, LP. Sigma brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.

Date: 05/2015  
Code: OFU